

Specification:

Standard Universal Screw



Size	P.N
Φ3.5×10	231481110C
Φ3.5×12	231481112C
Φ3.5×14	231481114C
Φ3.5×16	231481116C
Φ3.5×18	231481118C
Φ3.5×20	231481120C
Φ3.5×22	231481122C
Φ3.5×24	231481124C
Φ3.5×26	231481126C
Φ3.5×28	231481128C
Φ3.5×30	231481130C
Φ3.5×32	231481132C
Φ3.5×34	231481134C
Φ3.5×36	231481136C
Φ3.5×38	231481138C
Φ3.5×40	231481140C
Φ3.5×42	231481142C
Φ3.5×44	231481144C
Φ3.5×46	231481146C
Φ3.5×48	231481148C

Connecting Rod



Size	P.N
Φ3.5×40	231489040
Φ3.5×60	231489060
Φ3.5×80	231489080
Φ3.5×100	231489100
Φ3.5×120	231489120
Φ3.5×160	231489160
Φ3.5×200	231489200
Φ3.5×300	231489300
Φ3.5×400	231489400

Pre-bent Connecting rod



Size	P.N
Φ3.5: 85/120	231488120
Φ3.5: 85/140	231488140
Φ3.5: 85/160	231488160
Φ3.5: 85/200	231488200

Adjustable Connecting Rod



Size	P.N
Φ3.5 72/120	231489112
Φ3.5 72/220	231489122

Crosslink



Size	P.N
21	231480121C
28	231480128C
31	231480131C
35	231480135C
42	231480142C
49	231480149C
56	231480156C

Transconnector Cage Nuts



Size	P.N
/	231480004C

Transverse Connector Compression Screw



Size	P.N
/	231480003C

Occipital Plate



Size	P.N
27	231480027
31	231480031
37	231480037
45	231480045

Occipital Fixation Screw



Size	P.N
Φ4.5×6	231480706C
Φ4.5×8	231480708C
Φ4.5×10	231480710C
Φ4.5×12	231480712C

Outboard Shift Connector



Size	P.N
/	231480005

Axial Rod/ Rod Connector



Size	P.N
/	231480100

Occipital Plate Connector



Size	P.N
/	231480000

Compression Set Screw



Size	P.N
2.5	231480001C
3.2	231480002C

Titanium Cable Connector



Size	P.N
/	231480006

Lateral Rod/ rod connectors



Size	P.N
/	231480101

Leading the way to a healthy world

Canwell Medical Co.,Ltd

Add: No.466 South Xianhua Street
CN, Zhejiang Province, Jinhua City
321016

Web: <http://www.canwell.com.cn>



Spine

**CanPCS Occipitocervical
Posterior Fixation System**



CanPCS Occipitocervical Posterior Fixation System


Features and Advantages:

- Universal screw 60° variable Angle design, larger screw entry Angle
- Equipped with multi-diameter combined rods to meet different surgical needs of occipitocervical thoracic surgery
- A variety of implant options, complete specifications and models
- Special transverse connector makes operation more convenient
- Ultra small size, low profile excellent performance
- Color oxidation, with different colors to distinguish

Multiple Specifications of Occipital Plate:

27mm, 31mm, 37mm, 45mm occipital plate, with occipital connector


Occipital fixation screw:
length 6-16mm ϕ 4.5mm ϕ 5.25MM selection



Pedicle Screw:

Polyaxial screws: Standard Polyaxial, offset universal, Long shank universal

- Different types of screw color distinction
- Screw length: 8-50mm
- Screw diameter: 3.5mm, 4.0mm, 4.35mm



Multiple Forms of Titanium Rod

Standard connection rod, Pre-bent rod, Adjustable rod, Double diameter rod



Crosslink

Two options



Head-Head Cross Connector



Low profile Cross Connector

Complete Accessories

Can be seamlessly connected with the thoracolumbar posterior screw system

