


Various Cartridges



D-Series Cartridges

Model	Application	Staple height (mm)	Closed height(mm)	Staple quantity (pc)	Staple line(mm)	Color define
DCR46M	Septum	2.0	0.75	70	46	Gray
DCR46W	Thin tissue	2.5	1.0	70	46	White
DCR46B	Normal tissue	3.5	1.5	70	46	Blue
DCR46Y	Normal/thick tissue	3.8	1.8	70	46	Yellow
DCR46G	Thick tissue	4.1	2.0	70	46	Green
DCR60M	Septum	2.0	0.75	88	60	Gray
DCR60W	Thin tissue	2.5	1.0	88	60	White
DCR60B	Normal tissue	3.5	1.5	88	60	Blue
DCR60Y	Normal/thick tissue	3.8	1.8	88	60	Yellow
DCR60G	Thick tissue	4.1	2.0	88	60	Green
DCR60H	Extra Thick tissue	4.4	2.3	88	60	Black
D-Series Cartridges :				6 Rows		

T-Series Cartridges(with Non-slip technology)

Model	Application	Staple height (mm)	Closed height(mm)	Staple quantity (pc)	Staple line(mm)	Color define
TCR46M	Septum	2.0	0.75	70	46	Gray
TCR46W	Thin tissue	2.5	1.0	70	46	White
TCR46B	Normal tissue	3.5	1.5	70	46	Blue
TCR46Y	Normal/thick tissue	3.8	1.8	70	46	Yellow
TCR46G	Thick tissue	4.1	2.0	70	46	Green
TCR60M	Septum	2.0	0.75	88	60	Gray
TCR60W	Thin tissue	2.5	1.0	88	60	White
TCR60B	Normal tissue	3.5	1.5	88	60	Blue
TCR60Y	Normal/thick tissue	3.8	1.8	88	60	Yellow
TCR60G	Thick tissue	4.1	2.0	88	60	Green
TCR60H	Extra Thick tissue	4.4	2.3	88	60	Black
T-Series Cartridges :				6 Rows		



EnDrive™ Disposable Powered Endoscopic Linear Cutting Staplers



Ningbo Verykind Medical Device Co., Ltd

#100 Jinghua Road High-tech industrial development zone, Ningbo 315040 Zhejiang China


Contact us

Tel: +86-574-87802005 Fax: +86-574-87910671 www.nbverykind.com

Email: bettybian@nbverykind.com

Ningbo Verykind Medical Device Co., Ltd


Design of 60° articulation angle which can help the stapler reaches the pelvic easily with better solutions for low and ultra-low rectal cancer resection.



60° Articulation Angle

- Easily reaches the narrow space or even the ultra-low pelvic rectum resection vertically.
- Offers better visibility during laparoscopic surgery.
- Has big advantage of radical resection of low and ultra-low rectal carcinoma.
- Keeps benign tissues as much as possible during laparoscopic surgery.
- Reduces the occurrence of anastomotic fistula.
- Through 60° articulation to increase the possibilities of reserving anus operation in rectal cancer to improve patients' life quality.

PRODUCT ADVANTAGES



The design of 60° articulation angle, guarantee the operation of the low and ultra-low rectal cancer resection.


The innovative 6 rows stair-stepping, designed staple formation offers a superior haemostatic effect.

To prevent mistakenly firing.

5 Invention Patents 13 Utility Model Patents

Battery motivated stapling to guarantee smoothly firing and operation safety

Advanced automatically stapling function, ensuring smoothly and stably firing during the surgical procedure. Effectively prevent the tissue being pushed out of compressed cartridge by the firing force of manual endoscopic stapler. The knife is reserved automatically after firing in order to avoid stuck during pulling back. It's capable of offering average cut and suture as well as the excellent stapling formation.



Non-slip technology
(For T-Series)

Non-slip surface has better tissue compression for better anastomosis result.

Battery motivated

Bring more stable operation.

Item No.	Model	Shaft length (mm)
1	DCE46	280
2	DSE46	340
3	DLE46	440
4	DCE60	280
5	DSE60	340
6	DLE60	440